



hCG Rapid Pregnancy Test



Key Features:

- Simple one-step test
- Easy to use & read
- Accuracy: 99%
- Read results within 5 minutes
- Detects 25mIU/ml hCG
- Room temperature storage

Specifications:

- 1) Specimen: urine
- 2) Format: cassette
- 3) Sensitivity: 25mIU/mL
- 4) One kit includes 25 tests and pipettes in individual foil pouches
- 5) Pouches rinted with 'How To Read' Diagrams

Ordering Details:

- Cassettes: 25 tests
- Reference: RSFHC102-25

- ▶ The Suresign hCG Rapid Pregnancy test is a rapid test to detect the presence of hCG in urine specimens in a qualitative format sensitive to 25mIU hCG/ml.
- ▶ The test utilizes a combination of monoclonal and polyclonal antibody reagents to selectively detect elevated level of hCG in urine or serum.
- ▶ The immunological specificity of the test virtually eliminates cross reactivity interferences from the structurally related glycoprotein hormones hFSH, hLh, hTSH at physiological levels.

Award-winning suppliers to the NHS supply chain in Great Britain

CIGA Healthcare Ltd,
Kildowney House, Pennybridge Industrial Estate, Ballymena, BT42 3HB, UK
Tel: +44 (0) 2825 685 385
Fax: +44 (0) 2825 685 145
Email: info@cigahealthcare.com
Web: www.suresign.com

