



PULSE OXIMETER FINGERTIP TYPE

Measure: SpO₂ & PR



FEATURES OF THE JUMPER MEDICAL FINGERTIP PULSE OXIMETER

- **COLOUR OLED DISPLAY WITH 4 DISPLAY DIRECTIONS**
- **DETACHABLE LANYARD**
- **SpO₂**
 - MEASUREMENT RANGE: 35-100%
- **PULSE RATE:**
 - MEASUREMENT RANGE: 25-250BPM
 - PULSE INTENSITY: BARGRAPH INDICATOR
- **POWER REQUIREMENT: TWO AAA ALKALINE BATTERIES**
- **FEATURES:**
 - CE MARKED FOR CLINICAL & HOME USE
 - ACCURATELY MEASURES SpO₂ & PR
 - SpO₂ & PR DISPLAY, PR WAVEFORM & BARGRAPH DISPLAY
 - LOW BATTERY INDICATOR
 - FOUR DIRECTIONS & SIX DISPLAY MODES
 - AUTOMATICALLY POWERS OFF WHEN NOT IN USE
 - PULSE SOUND FUNCTION



CIGA Healthcare Ltd
Kildowney House, Pennybridge Industrial Estate, Ballymena BT42 3HB

suresignprofessional.com | cigahealthcare.com

Tel: 028 2568 5385 Email: info@cigahealthcare.com
SB-424-548A