



SURESIGN FINECARE POINT OF CARE D-DIMER TEST

Aids in the management of DVT & PE



The Suresign Finecare D-dimer is a rapid in vitro diagnostic test for the quantitative determination of D-Dimer in human blood. The point of care D-Dimer test helps guide decisions in the diagnosis of urgent situations such as DVT and PE. The Suresign Finecare D-dimer test fits well with the requirements of D-Dimer testing in the NICE Guideline CG144

- **RAPID POINT OF CARE TEST – RESULTS IN 5 MINUTES**
- **SIMPLE TO USE TWO-STEP TEST**
- **REQUIRES A SMALL VOLUME OF VENOUS BLOOD 15µL OR PLASMA 18µL**
- **MEASURING RANGE 0.1-10mg/L**
- **25 TEST KIT WITH LONG SHELF LIFE**
- **ROOM TEMPERATURE STORAGE**
- **FLUORESCENT SIGNALLING ALLOWING HIGH SENSITIVITY AND SPECIFICITY**
- **PROVIDES LABORATORY STANDARD RESULTS**
- **NO CALIBRATION REQUIRED**



CIGA Healthcare Ltd
Kildowney House, Pennybridge Industrial Estate, Ballymena BT42 3HB

suresignprofessional.com | cigahealthcare.com

Tel: 028 2568 5385 Email: info@cigahealthcare.com